

Product datasheet for PH315131

ANKRD16 (NM_001009941) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	ANKRD16 MS Standard C13 and N15-labeled recombinant protein (NP_001009941)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC215131
Predicted MW:	39.3 kDa
Protein Sequence:	>RC215131 protein sequence Red=Cloning site Green=Tags(s)
	<p>MAQPGDPRRLCRLVQEGRLRALKEELQAAGGCPGPAGDTLLHCAARHGHRDVLAYLAEAWGMDIEATNRD YKRPLHEAASMGRDCVRYLLGRGAADVCLKKADWTPLMMACTRKNLGVIQELVEHGANPLLKNKDGWNS FHIASREGDPLILQYLLTVCPGAWKTESKIRRTPLHTAAMHGHEAVKVLLKRCQYEPDYRDNCGVTALM DAIQCGHIDVARLLLDEHGACLSAEDSLGAQALHRAAVTGQDEAIRFLVSELGVDVVDVRATSTHLTALHY AAKEGHTSTIQTLTLLSLGADINSKDEKNRSALHLACAGQHLACAKFLLQSGLKDSEDITGTLAQQLPRRAD VLQGS GHSAMT</p> <p>TRTRPLEQKLI SEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_001009941</u>
RefSeq Size:	1683
RefSeq ORF:	1083
Locus ID:	54522
UniProt ID:	<u>Q6P6B7</u>




[View online »](#)

Cytogenetics: 10p15.1

Summary: Required to prevent the misactivation of serine (Ser) with tRNA(Ala) by promoting the hydrolysis of Ser-mischarged tRNA(Ala), thereby playing a role in translational fidelity. Binds directly to the catalytic domain of AARS/AlaRS and captures Ser that is misactivated by AARS/AlaRS, preventing the charging of Ser adenylates to tRNA(Ala) and precluding Ser misincorporation in nascent peptides.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified ANKRD16 protein (Cat# [TP315131]). The protein was produced from HEK293T cells transfected with ANKRD16 cDNA clone (Cat# [RC215131]) using MegaTran 2.0 (Cat# [TT210002]).