

Product datasheet for PH315180

NFAT1 (NFATC2) (NM_012340) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	NFATC2 MS Standard C13 and N15-labeled recombinant protein (NP_036472)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC215180
Predicted MW:	99.6 kDa
Protein Sequence:	>RC215180 representing NM_012340 Red=Cloning site Green=Tags(s)

MNAPERQPQPDGGDAPGHEPGGSPQDELDFSILFDYEYLNPNEEEEPNAHKVASPPSGPAYPDDVLDYGLK
PYSPLASLSGEPGRFGEPDRVGPQKFLSAAKPAGASGLSPRIEITPSHELIQAVGPLRMRDAGLLVEQP
PLAGVAASPRFTLPVPGFEGYREPLCLSPASSGSSASFISDTFSPYTPSPCVSPNNGPDDLCQFQNI
HYSPTSPIMSPRTSLAEDSCLGRHSPVPRPASRSSSPGAKRRHSCEALVALPPGASPQRSRSPSPQPS
SHVAPQDHGSPAGYPPVAGSAVIMDALNSLATDSPCGIPPKMWTSPDPSVSAAPSKAGLPRHIYPAVE
FLGPCEQGERRNSAPESILLVPPTWPKPLVPAIPICSIPTASLPLEWPLSSQSGSYELRIEVQPKPHH
RAHYETEGSRGAVKAPTGGHPVVQLHGYMENKPLGLQIFIGTADERILKPHAFYQVHRITGKTVTTTSYE
KIVGNTKVL EIPLEPKNNMRATIDCAGILKLRNADIELRKGETDIGRKNTRVRLVFRVHIPESSGRIVSL
QTASNPIECSQRSHELPMVERQDTSCLVYGGQMILTGQNF TSESKVVFTEKTTDGGQIWEMEATVDK
DKSQPNMLFVEIPEYRNKHIRTVPKVN FYINGKRKRSQPQHFTYHPVPAIKTEPTDEYDPTLICSPTHG
GLGSQPYYPQHPMVAESPSCLVATMAPCQFRTGLSSPDARYQQNPAAVLYQRSKSLSPSLLGYQQPAL
MAAPLSLADAHRVSVL VHAGSQGQSSALLHPSPTNQQASPV IHYSP TNQQLRGSHQEFQHIMYCENFAPG
TTRPGPPPVSQGRQLSPGSYPTVIQQNATSQRAAKNGPPVSDQKEVLPAGVTIKQEQLNDQTYLDEL I
DTHLSWIQNIL

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.



[View online »](#)

RefSeq: [NP_036472](#)

RefSeq Size: 3254

RefSeq ORF: 2763

Synonyms: NFAT1; NFATP

Locus ID: 4773

UniProt ID: [Q13469](#)


Cytogenetics: 20q13.2

Summary: This gene is a member of the nuclear factor of activated T cells (NFAT) family. The product of this gene is a DNA-binding protein with a REL-homology region (RHR) and an NFAT-homology region (NHR). This protein is present in the cytosol and only translocates to the nucleus upon T cell receptor (TCR) stimulation, where it becomes a member of the nuclear factors of activated T cells transcription complex. This complex plays a central role in inducing gene transcription during the immune response. Alternate transcriptional splice variants encoding different isoforms have been characterized. [provided by RefSeq, Apr 2012]

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Axon guidance, B cell receptor signaling pathway, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, T cell receptor signaling pathway, VEGF signaling pathway, Wnt signaling pathway

Product images:



Coomassie blue staining of purified NFATC2 protein (Cat# [TP315180]). The protein was produced from HEK293T cells transfected with NFATC2 cDNA clone (Cat# [RC215180]) using MegaTran 2.0 (Cat# [TT210002]).