

Product datasheet for PH315613

TATDN3 (NM_001042553) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	TATDN3 MS Standard C13 and N15-labeled recombinant protein (NP_001036018)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC215613
Predicted MW:	30.3 kDa
Protein Sequence:	>RC215613 protein sequence Red=Cloning site Green=Tags(s) MRAAGVGLVDCHCHLSAPDFDRDLDDVLEKAKKANVVALVAVAEHSGEFEKIMQLSERYNGFVLPCLGVH PVQGLPPEDQRSVTLKDLVALPIIENYKDRLLAIGEVGLDFSPRFAGTGEQKEEQRVLIRQIQLAKRL NLPVNVHSRSAGRPTINLLQEQGAEKVLLHAFDGRPSVAMEGVRAGYFFSIPPSIIRSGQKQKLVKQLPL TSICLETDSPALGPEKQVRNEPWNISISAEYIAQVKGISVEEVIEVTTQNALKLFPKLRHLLQK TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001036018
RefSeq Size:	2585
RefSeq ORF:	822
Locus ID:	128387
UniProt ID:	Q17R31
Cytogenetics:	1q32.3



[View online »](#)

Summary: Putative deoxyribonuclease.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified TATDN3 protein (Cat# [TP315613]). The protein was produced from HEK293T cells transfected with TATDN3 cDNA clone (Cat# [RC215613]) using MegaTran 2.0 (Cat# [TT210002]).