

Product datasheet for PH315846

MAGIX (NM_001099682) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	MAGIX MS Standard C13 and N15-labeled recombinant protein (NP_001093152)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC215846
Predicted MW:	26.5 kDa
Protein Sequence:	>RC215846 representing NM_001099682 Red =Cloning site Green =Tags(s) MEPRTGGAANPKGSRGSRGPSPLAGPSARQLLARLDARPLAARAAYDVAALVRRAGATLRLRKEAVSVL DSADIEVTDSRLPHATIVDHRPQVGDVLHINGESTQGLTHAQAVIRIRAGGPQLHLVIRPLETHPGKP RGVGEPRKGVDRSPDPGGPEVTGSRSSSTSLVQHPPSRTTLKKTRGSEPSPEAAADGPTVSPERRAED PNDQIPGSPGPWLVPSEERLSRALGVRGAAQLAQEMAAGRHH TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001093152
RefSeq Size:	1907
RefSeq ORF:	759
Synonyms:	JM10; PDZX
Locus ID:	79917
UniProt ID:	Q9H6Y5 , A0A087X263



[View online »](#)

Cytogenetics: Xp11.23

Product images:



Coomassie blue staining of purified MAGIX protein (Cat# [TP315846]). The protein was produced from HEK293T cells transfected with MAGIX cDNA clone (Cat# [RC215846]) using MegaTran 2.0 (Cat# [TT210002]).