

Product datasheet for PH316875

DPP6 (NM_001936) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	DPP6 MS Standard C13 and N15-labeled recombinant protein (NP_001927)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC216875
Predicted MW:	91.2 kDa
Protein Sequence:	>RC216875 representing NM_001936 Red=Cloning site Green=Tags(s)

MTTAKEPSASGKSVQQEQELVGSNPPQRNWKGIALLVILVICSLIVTSVILLTPAEDNSLSQKKKVT
VEDLFSEDFKIHDPKAWISDTEFIYREQKGTVRLWVETNTSTVLEGGKIESLRAIRYIEISPDREYAL
FSYNVEPIYQHSYTGYYVLSKIPHGDPQSLDPPEVSNAKLQYAGWGPKGQQLIFIFENNIYYCAHVKGQA
IRVVSTGKEGVIYNGLSDWLYEEEILKTHIAHWSPDGTRELAYAAINDSRVPIMELPTYTGSYPTVKPY
HYPKAGSENPISLHVIGLNGPTHLEMMPPDDPRMREYYITMVKATSTKVAVTWNRAQNVSILTLCD
ATTGVCTKKHEDESEAWLHRQNEEPVFSKDKRKFIRAIPOGGRGKFYHITVSSSQPNSSNDNIQSITS
GDWDVTKILAYDEKGNKIYFLSTEDLPRRRQLYSANTVGNFNRQCLSCDLVENCTYFSASFSSMDFFLL
KCEGPGVPMVTVHNTTDDKKMFDETNEHVKKAINDRQMPKVEYRDIEIDYDNLPMQILKPATFTDTHY
PLLLVVDGTPGSQVAEKFEVSWETVMVSSHGAVVVKCDGRSGFGTKLLHEVRRRLGLLEEKDQMEAV
RMLKEQYIDRTRVAVFGKDYGGYLSTYILPAKGENQGGTFTCGSALSPITDFKLYASAFSERYLGLHGL
DNRAYEMTKVAHRVSALEEQQFLIIHPTADEKIHQHTAELITQLIRGKANYSLQIYPDESHYFTSSSLK
QHLYRSIINFFVECFRIQDKLLTVTAKEDDEED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001927




[View online >](#)

RefSeq Size:	3560
RefSeq ORF:	2409
Synonyms:	DPL1; DPPX; MRD33; VF2
Locus ID:	1804
UniProt ID:	A7E2E4
Cytogenetics:	7q36.2

Summary: This gene encodes a single-pass type II membrane protein that is a member of the peptidase S9B family of serine proteases. This protein has no detectable protease activity, most likely due to the absence of the conserved serine residue normally present in the catalytic domain of serine proteases. However, it does bind specific voltage-gated potassium channels and alters their expression and biophysical properties. Variations in this gene may be associated with susceptibility to amyotrophic lateral sclerosis and with idiopathic ventricular fibrillation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014]

Protein Families: Druggable Genome, Protease, Transmembrane

Product images:



Coomassie blue staining of purified DPP6 protein (Cat# [TP316875]). The protein was produced from HEK293T cells transfected with DPP6 cDNA clone (Cat# [RC216875]) using MegaTran 2.0 (Cat# [TT210002]).