

Product datasheet for PH316883

G protein alpha S (GNAS) (NM_001077489) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	GNAS MS Standard C13 and N15-labeled recombinant protein (NP_001070957)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC216883
Predicted MW:	44 kDa
Protein Sequence:	>RC216883 representing NM_001077489 Red =Cloning site Green =Tags(s)

MGCLGNSKTEDQRNEEKAQREANKKIEKQLQKDKQVYRATHRLLLGAGESGKSTIVKQMRILHVNFGNG
 DEKATKVQDIKNNLKEAIETIVAAMSNLVPPVELANPENQFRVDYILSVMNVPDFDPPEFYEHAKALWE
 DEGVRACYERSNEYQLIDCAQYFLDKIDVIKQADYVPSDQDLLRCRVLTSIGIFETKFQVDKVNFMFDVG
 GQRDERRKWIQCNDVTAIIFVVAASSYNMVIREDNQTNRLQEALNLFKSIWNNRWLRTISVILFLNKQD
 LLAEKVLGAKSKIEDYFPEFARYTTPEDATPEPGEDPRVTRAKYFIRDEFRLISTASGDGRHYCYPHFTC
 AVDTENIRRVFNDCRDIIQRMHLRQYELL

TRRLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_001070957</u>
RefSeq Size:	1881
RefSeq ORF:	1137
Synonyms:	AHO; C20orf45; GNAS1; GPSA; GSA; GSP; NESP; PITA3; POH; SCG6; SgVI
Locus ID:	2778


[View online »](#)

UniProt ID: [Q5FWY2](#), [O95467](#), [P63092](#)


Cytogenetics: 20q13.32

Summary: This locus has a highly complex imprinted expression pattern. It gives rise to maternally, paternally, and biallelically expressed transcripts that are derived from four alternative promoters and 5' exons. Some transcripts contain a differentially methylated region (DMR) at their 5' exons, and this DMR is commonly found in imprinted genes and correlates with transcript expression. An antisense transcript is produced from an overlapping locus on the opposite strand. One of the transcripts produced from this locus, and the antisense transcript, are paternally expressed noncoding RNAs, and may regulate imprinting in this region. In addition, one of the transcripts contains a second overlapping ORF, which encodes a structurally unrelated protein - Alex. Alternative splicing of downstream exons is also observed, which results in different forms of the stimulatory G-protein alpha subunit, a key element of the classical signal transduction pathway linking receptor-ligand interactions with the activation of adenylyl cyclase and a variety of cellular responses. Multiple transcript variants encoding different isoforms have been found for this gene. Mutations in this gene result in pseudohypoparathyroidism type 1a, pseudohypoparathyroidism type 1b, Albright hereditary osteodystrophy, pseudopseudohypoparathyroidism, McCune-Albright syndrome, progressive osseous heteroplasia, polyostotic fibrous dysplasia of bone, and some pituitary tumors. [provided by RefSeq, Aug 2012]

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Calcium signaling pathway, Dilated cardiomyopathy, Gap junction, GnRH signaling pathway, Long-term depression, Melanogenesis, Taste transduction, Vascular smooth muscle contraction, Vibrio cholerae infection

Product images:



Coomassie blue staining of purified GNAS protein (Cat# [TP316883]). The protein was produced from HEK293T cells transfected with GNAS cDNA clone (Cat# [RC216883]) using MegaTran 2.0 (Cat# [TT210002]).