

Product datasheet for PH317380

FAM156B (NM_001099684) Human Mass Spec Standard


Product data:

Product Type:	Mass Spec Standards
Description:	FAM156B MS Standard C13 and N15-labeled recombinant protein (NP_001093154)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC217380
Predicted MW:	24.2 kDa
Protein Sequence:	>RC217380 representing NM_001099684 Red=Cloning site Green=Tags(s) MDPLQKRNPASPSKSSPMTAAETSQEGPAPSQPSYSEQPMMGLSNLSPGPGHSQAVPLPEGLLHQRYREE KTLERRWERLEFLQRKKAFLRHVRRRHRDHMAPYAVGREARISPLGDRSQNRFCQCECRYCQSHRPNLSG IPGESNRAPHSSWETLVQGLSGLTSLGTNQPGPLPEAALQPQETEEKRQERQQESKIMFQRLKQWL EEN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001093154
RefSeq Size:	1701
RefSeq ORF:	639
Synonyms:	TMEM29B
Locus ID:	727866
UniProt ID:	Q8NDB6


[View online »](#)

Cytogenetics: Xp11.22

Product images:



Coomassie blue staining of purified FAM156B protein (Cat# [TP317380]). The protein was produced from HEK293T cells transfected with FAM156B cDNA clone (Cat# [RC217380]) using MegaTran 2.0 (Cat# [TT210002]).