

Product datasheet for PH317667

LDLRAD2 (NM_001013693) Human Mass Spec Standard


Product data:

Product Type:	Mass Spec Standards
Description:	LDLRAD2 MS Standard C13 and N15-labeled recombinant protein (NP_001013715)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC217667
Predicted MW:	28.4 kDa
Protein Sequence:	>RC217667 representing NM_001013693 Red=Cloning site Green=Tags(s) MEACLLQLPQRLLLLGAAALTATALETADLAELCGQTWQDGLLLRSHAASRRFYFVAPDTCGLWVQA AAPGDRIRFQFRFFLVYSLTPAPPALNTSSPAPADPCAPGSYLRFYEGPPGAPRPLGSPLCGLNIPVPVA SSGPFLGLRLVTRGRQPRVDFVGEVTSFRLGPCGAYFRCQNGRCIPSSLVCDPWGMDNCGDQSGSWSP ADCRGSPVPSQTGSTAHTSRSLTPSPALGSAGSLWIAAERSSPAGRDPTRQDAALEGSTE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001013715
RefSeq Size:	4012
RefSeq ORF:	816
Locus ID:	401944
UniProt ID:	Q5SZ11
Cytogenetics:	1p36.12


[View online »](#)

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified LDLRAD2 protein (Cat# [TP317667]). The protein was produced from HEK293T cells transfected with LDLRAD2 cDNA clone (Cat# [RC217667]) using MegaTran 2.0 (Cat# [TT210002]).