

Product datasheet for PH318920

Glucokinase (GCK) (NM_033508) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	GCK MS Standard C13 and N15-labeled recombinant protein (NP_277043)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC218920
Predicted MW:	51.9 kDa
Protein Sequence:	>RC218920 representing NM_033508 Red=Cloning site Green=Tags(s) MPRPRSQLPQPNSQVEQILAEFQLQEEDLKKVMRRMQKEMDRGLRLETHEEASVKMLPTYVRSTPEGSEV GDFLSLDLGGTNFRVMLVKVGEEGEQWSVKTKHQMYSIPEDAMTGAEMLFDYISECISDFLDKHQMKH KKLPLGFTFSFPVRHEDIDKGILLNWTGKFKASGAEGNNVVGLLRDAIKRRGDFEMDVVAMVNDTVATMI SCYYEDHQCEVGMIVGTGCNACYMEEMQNVELVEGDEGRMCVNTIEWGAFGDSGELDEFLLDYDRLVDESS ANPGQQLYEKLIIGGKYMVELVRLVLLRLVDENLLFHGEASEQLRTRGAFETRFVSVQVESDTGDRKQIYNI LSTLGLRPSTTDCDIVRRACESVSTRAAHMCSAGLAGVINRMRESRSEDVMRITVGVGDGSVYKLLHPSFKE RFHASVRRRLTPSCEITFIESEEGSGRGAALVSAVACKKACMLGQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_277043</u>
RefSeq Size:	2566
RefSeq ORF:	1392
Synonyms:	FGQTL3; GK; GLK; HHF3; HK4; HKIV; HXKP; LGLK; MODY2; PNDM1



[View online »](#)

Locus ID: 2645

UniProt ID: [P35557](#), [P35557-3](#)

Cytogenetics: 7p13

Summary: This gene encodes a member of the hexokinase family of proteins. Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. In contrast to other forms of hexokinase, this enzyme is not inhibited by its product glucose-6-phosphate but remains active while glucose is abundant. The use of multiple promoters and alternative splicing of this gene result in distinct protein isoforms that exhibit tissue-specific expression in the pancreas and liver. In the pancreas, this enzyme plays a role in glucose-stimulated insulin secretion, while in the liver, this enzyme is important in glucose uptake and conversion to glycogen. Mutations in this gene that alter enzyme activity have been associated with multiple types of diabetes and hyperinsulinemic hypoglycemia. [provided by RefSeq, Aug 2017]

Protein Families: Druggable Genome

Protein Pathways: Amino sugar and nucleotide sugar metabolism, Galactose metabolism, Glycolysis / Gluconeogenesis, Insulin signaling pathway, Maturity onset diabetes of the young, Metabolic pathways, Starch and sucrose metabolism, Type II diabetes mellitus

Product images:



Coomassie blue staining of purified GCK protein (Cat# [TP318920]). The protein was produced from HEK293T cells transfected with GCK cDNA clone (Cat# [RC218920]) using MegaTran 2.0 (Cat# [TT210002]).