

Product datasheet for PH318995

TPM1 (NM_001018020) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	TPM1 MS Standard C13 and N15-labeled recombinant protein (NP_001018020)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC218995
Predicted MW:	32.6 kDa
Protein Sequence:	>RC218995 representing NM_001018020 Red=Cloning site Green=Tags(s) MDAIIKKKMQMLKLDKENALDRAEQAEADKKAEDRSKQLEEDIAAKEKLLRVSEDERDRVLEELHKAEDS LLAAEEAAAKAEADVASLNRRISLVEEELDRAQERLATALQKLEEAKEAADESERGMKVIESRAQKDEEK MEIQEIQLKEAKHIAEDADRKYEEVARKLVIIIESDLERAEEAELSEGQVRQLEEQRLIMDQTLKALMAA EDKYSQKEDRYEEEIKVLSDKLKEAETRAEFAERSVTKLEKSIDDLEEKVAHAKEENLSMHQMLDQTLLE LNNM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001018020
RefSeq Size:	1797
RefSeq ORF:	852
Synonyms:	C15orf13; CMD1Y; CMH3; HEL-S-265; HTM-alpha; LVNC9; TMSA
Locus ID:	7168



[View online »](#)

UniProt ID: [P09493](#), [Q15513](#)


Cytogenetics: 15q22.2

Summary: This gene is a member of the tropomyosin family of highly conserved, widely distributed actin-binding proteins involved in the contractile system of striated and smooth muscles and the cytoskeleton of non-muscle cells. Tropomyosin is composed of two alpha-helical chains arranged as a coiled-coil. It is polymerized end to end along the two grooves of actin filaments and provides stability to the filaments. The encoded protein is one type of alpha helical chain that forms the predominant tropomyosin of striated muscle, where it also functions in association with the troponin complex to regulate the calcium-dependent interaction of actin and myosin during muscle contraction. In smooth muscle and non-muscle cells, alternatively spliced transcript variants encoding a range of isoforms have been described. Mutations in this gene are associated with type 3 familial hypertrophic cardiomyopathy. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM)

Product images:



Coomassie blue staining of purified TPM1 protein (Cat# [TP318995]). The protein was produced from HEK293T cells transfected with TPM1 cDNA clone (Cat# [RC218995]) using MegaTran 2.0 (Cat# [TT210002]).