

Product datasheet for PH319023

hnRNP K (HNRNPK) (NM_031263) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	HNRNPK MS Standard C13 and N15-labeled recombinant protein (NP_112553)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC219023
Predicted MW:	51 kDa
Protein Sequence:	>RC219023 protein sequence Red=Cloning site Green=Tags(s)

METEQPEETFPNTEETNGEFGKRPAMEDMEEEQAFKRSRNTDEMVELRILLQSKNAGAVIGKGGKNIKALRT
DYNASVSVPDSSGPERILSISADIETIGEILKKIIPBLEEGLQLPSPTATSQLPLESDAVECLNYQHYKG
SDFDCELRLLIHQSLAGGIIGVKGAKIKELRENTQTTIKLFQECCPHSTDRVVLIGGKPDVVECIKIIIL
DLISESPIKGRAQPYDPNFYDETYDYGFTMMFDDRRGRPVGFPMRGRGGFDRMPPGRGGRPMPSSRRDY
DDMSPRRGPPPPPPGRGGSRARNLPLPPPPPRGGDLMA YDRRGRPGDRYDGMVGFSADETWDS AID
TWSPSEWQMAYEPQGGSGYDYSYAGGRGSYDLGGPIITQVTIPKDLAGSIIGKGGQRIKQIRHESGAS
IKIDEPLEGSEDRIITITGTQDQIQNAQYLLQNSVKQYADVEGF

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<u>NP_112553</u>
RefSeq Size:	2960
RefSeq ORF:	1392
Synonyms:	AUKS; CSBP; HNRPK; TUNP




[View online »](#)

Locus ID: 3190
UniProt ID: [P61978](#)
Cytogenetics: 9q21.32

Summary: This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene is located in the nucleoplasm and has three repeats of KH domains that binds to RNAs. It is distinct among other hnRNP proteins in its binding preference; it binds tenaciously to poly(C). This protein is also thought to have a role during cell cycle progression. Several alternatively spliced transcript variants have been described for this gene, however, not all of them are fully characterized. [provided by RefSeq, Jul 2008]

Protein Pathways: Spliceosome

Product images:



Coomassie blue staining of purified HNRNPK protein (Cat# [TP319023]). The protein was produced from HEK293T cells transfected with HNRNPK cDNA clone (Cat# [RC219023]) using MegaTran 2.0 (Cat# [TT210002]).