

Product datasheet for PH319629

CYLD (NM_015247) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	CYLD MS Standard C13 and N15-labeled recombinant protein (NP_056062)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC219629
Predicted MW:	107.1 kDa
Protein Sequence:	>RC219629 representing NM_015247 Red=Cloning site Green=Tags(s)

MSSGLWSQEKVTSFYWEERIFYLLQLQECVTDKQTQKLLKVPKGSIGQYIQDRSVGHSRIPSAKGGKKNQI
GLKILEQPHAVLFDVDEKDVVEINEKFTLELLAITNCEERFSLFKNRNRLSKGLQIDVGCVPKVQLRSGEE
KFPGVVFRFRGPLLAERTVSGIFFVELLEEGRGQGFDTGVYQGGKQLFQCDEDCGVFVALDKLELIEDDDT
ALESDYAGPGDTMQVELPPLEINSRVSLKVGETIESGTVIFCDVLPGKESLGYFVGVDMDNPIGNWDGRF
DGVQLCSFACVESTILLHINDIIPALSESVTQERRPPKLAFMSRGGVGDGSSSHNKPKATGSTSDPGNRN
RSELYFTLNGSSVDSQPQSKSKNTWYIDEVAEDPAKSLTEISTDFDRSSPPLQPPPVNSLTENRFHSLP
FSLTKMPNTNGSIGHSPLSLSAQSYMEELENTAPVQESPLAMPNGSHGLEVGSLAEVKENPPFYGVIRW
IGQPPGLNEVLAGLELEDEECAGCTDGTFRGTRYFTCALKKALFVKLKSRRPDSRFASLQPVSNQIERCNS
LAFGGYLSEVVEENTPPKMEKEGLEIMIGKKKGIQGHYNSCYLDSTLFCLFAFSSVLDTVLLRPKEKNDV
EYYSETQELLRTEIVNPLRIYGVYCATKIMKLRKILEKVEAASGFTSEEKDPPEFLNILFHILRVPELL
KIRSAGQKVQDCYFYQIFMEKNEKVGVPYTIQQLLEWSFINSNLKFAEAPSCLIIQMPRFKDFKFKKIF
PSLELNIIDLLEDTPRQCRI CGGLAMYECRECYDDPDISAGKIKQFCKTCNTQVHLHPKRLNHNKYNPVSL
PKDLPDWDWRHGCIPCQNMELFAVLCIETSHYVAFVKYKGDSSAWLFFDSMADRDRGGQNGFNIPQVTPCP
EVGEYLKMSLEDLHSLDSRRIQGCARRLLCDAYMCMYQSPTMSLYK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.



[View online »](#)


RefSeq:	NP_056062
RefSeq Size:	5371
RefSeq ORF:	2868
Synonyms:	BRSS; CDMT; CYLD1; CYLDI; EAC; FTDALS8; MFT; MFT1; SBS; TEM; USPL2
Locus ID:	1540
UniProt ID:	Q9NQC7
Cytogenetics:	16q12.1

Summary: This gene encodes a cytoplasmic protein with three cytoskeletal-associated protein-glycine-conserved (CAP-GLY) domains that functions as a deubiquitinating enzyme. Mutations in this gene have been associated with cylindromatosis, multiple familial trichoepithelioma, and Brooke-Spiegler syndrome. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: RIG-I-like receptor signaling pathway

Product images:



Coomassie blue staining of purified CYLD protein (Cat# [TP319629]). The protein was produced from HEK293T cells transfected with CYLD cDNA clone (Cat# [RC219629]) using MegaTran 2.0 (Cat# [TT210002]).