

Product datasheet for PH320343

MAX (NM_145113) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	MAX MS Standard C13 and N15-labeled recombinant protein (NP_660088)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC220343
Predicted MW:	11.9 kDa
Protein Sequence:	>RC220343 representing NM_145113 Red=Cloning site Green=Tags(s) MSDNDDEEVESDEEQPRFQSAADKRAHHNALERKRRDHDKDSFHSLRDSVPSLQGEKPQASRAQILDKAT EYIQYMRRNHHTHQQDIDDLKRQNALLEQQGESES TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_660088
RefSeq Size:	2169
RefSeq ORF:	315
Synonyms:	bHLHd4
Locus ID:	4149
UniProt ID:	P61244 , A0A024R682 , Q8TAX8
Cytogenetics:	14q23.3



[View online »](#)

Summary:


The protein encoded by this gene is a member of the basic helix-loop-helix leucine zipper (bHLHZ) family of transcription factors. It is able to form homodimers and heterodimers with other family members, which include Mad, Mxi1 and Myc. Myc is an oncoprotein implicated in cell proliferation, differentiation and apoptosis. The homodimers and heterodimers compete for a common DNA target site (the E box) and rearrangement among these dimer forms provides a complex system of transcriptional regulation. Mutations of this gene have been reported to be associated with hereditary pheochromocytoma. A pseudogene of this gene is located on the long arm of chromosome 7. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]

Protein Families:

Druggable Genome, Transcription Factors

Protein Pathways:

MAPK signaling pathway, Pathways in cancer, Small cell lung cancer

Product images:

Coomassie blue staining of purified MAX protein (Cat# [TP320343]). The protein was produced from HEK293T cells transfected with MAX cDNA clone (Cat# [RC220343]) using MegaTran 2.0 (Cat# [TT210002]).