

Product datasheet for PH320445

Liver Carboxylesterase 1 (CES1) (NM_001025195) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	CES1 MS Standard C13 and N15-labeled recombinant protein (NP_001020366)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC220445
Predicted MW:	62.59 kDa
Protein Sequence:	>RC220445 representing NM_001025195 Red=Cloning site Green=Tags(s)

MWLRAFILATLSASAAGHPSSPPVVDTVHGKVLGKFVSLEGFAQPVAIFLGIPFAKPPGLRFTPPQ
PAEPWSFVKNATSYPPMCTQDPKAGQLLSELFTRKENIPLKLSLSECLYLNIYTPADLTKKNRPLVMVWI
HGGGLMVGAASTYDGLALAAHENVVVVTIQYRLGIWGFSTGDEHSRGNWGHLDQVAALRWQDNIAISFG
GNPGSVTIFGESAGGESVSVLVSPLAKNLFHRAISESGVALTSVLVKKGDVKPLAEQIAITAGCKTTTS
AVMVHCLRQKTEEELLETTLKMKFLSLDLQGDPPRESQPLLGTVIDGMLLLKTPPELQAERNFHTVPYMG
INKQEFGLWIPQLMSYPLSEGQLDQKTAMSLWKSYPVLCIAKELIPEATEKYLGGTDDTVKKKDLFLD
LIADVMFGVPSVIVARNHRDAGAPTYMYEFQYRPSFSSDMKPKTVIGDHGDELFSVFGAPFLKEGASEEE
IRLSKMYMKFWANFARNGNPNGEGLPHWPEYNQKEGYLQIGANTQAAQKLDKEVAFWTNLFKAKAVEKP
PQTEHIEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001020366
RefSeq Size:	2027
RefSeq ORF:	1704



[View online »](#)

Synonyms: ACAT; CE-1; CEH; CES2; hCE-1; HMSE; HMSE1; PCE-1; REH; SES1; TGH

Locus ID: 1066

UniProt ID: [P23141](#)


Cytogenetics: 16q12.2

Summary: This gene encodes a member of the carboxylesterase large family. The family members are responsible for the hydrolysis or transesterification of various xenobiotics, such as cocaine and heroin, and endogenous substrates with ester, thioester, or amide bonds. They may participate in fatty acyl and cholesterol ester metabolism, and may play a role in the blood-brain barrier system. This enzyme is the major liver enzyme and functions in liver drug clearance. Mutations of this gene cause carboxylesterase 1 deficiency. Three transcript variants encoding three different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]

Protein Families: Druggable Genome

Protein Pathways: Drug metabolism - other enzymes

Product images:



Coomassie blue staining of purified CES1 protein (Cat# [TP320445]). The protein was produced from HEK293T cells transfected with CES1 cDNA clone (Cat# [RC220445]) using MegaTran 2.0 (Cat# [TT210002]).