


Locus ID: 91522

UniProt ID: [Q86Y22](#)

Cytogenetics: 5q35.3

Summary: COL23A1 is a member of the transmembrane collagens, a subfamily of the nonfibrillar collagens that contain a single pass hydrophobic transmembrane domain (Banyard et al., 2003 [PubMed 12644459]).[supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified COL23A1 protein (Cat# [TP320535]). The protein was produced from HEK293T cells transfected with COL23A1 cDNA clone (Cat# [RC220535]) using MegaTran 2.0 (Cat# [TT210002]).