

## Product datasheet for PH320559

### BARHL1 (NM\_020064) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	BARHL1 MS Standard C13 and N15-labeled recombinant protein (NP_064448)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC220559
Predicted MW:	34.9 kDa
Protein Sequence:	>RC220559 representing NM_020064 Red=Cloning site Green=Tags(s)  MEGSNGFGIDSILSHRAGSPALPKGDPLLGDCRSPLELSPRSESSDCSSPASPGRDCLETGTTPRGGAS GPGLDSHLQPGQLSAPAQSRVTSSFLIRDILADCKPLAACAPYSSSQPAAPEPGGRLAAKAAEDFRDK LDKSGSNASSDSEYKVEEGDREISSSRDPPVRLKKPRKARTAFTDHQLAQLERSFERQKYL SVQDRME LAASLNL TDTQVKTWYQNRRTKWKRQTAVGLELLAEAGNYSALQRMFSPYFYPQSLVSNLDPGAALYLY RGPSAPPALQRPLVPRI LIHGLQGASEPPPPLPPLAGVLPRAAQPR  TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_064448</a>
RefSeq Size:	1907
RefSeq ORF:	981
Locus ID:	56751
UniProt ID:	<a href="#">Q9BZE3</a>




[View online »](#)

Cytogenetics: 9q34.13

Protein Families: Transcription Factors

### Product images:



Coomassie blue staining of purified BARHL1 protein (Cat# [TP320559]). The protein was produced from HEK293T cells transfected with BARHL1 cDNA clone (Cat# [RC220559]) using MegaTran 2.0 (Cat# [TT210002]).