

Product datasheet for PH320978

ZNF483 (NM_001007169) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	ZNF483 MS Standard C13 and N15-labeled recombinant protein (NP_001007170)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC220978
Predicted MW:	29.4 kDa
Protein Sequence:	>RC220978 representing NM_001007169 Red=Cloning site Green=Tags(s) MQAVVPLNKMTAISPEPQTLASTEQNEVPRVVTSGEQEAILRGNAADAESFRQRFRWFCYSEVAGPRKAL SQLWELCNQWLRPDIHTKEQILELLVFEQFLTILPGEIRIWKVKSQHPESSEEVVTLIEDLTQMLEEKDPV SQDSTVSQEENSKEDKMVTVCPNTESCESITLKDVAVNFSRGEWKLEPFQKELYKEVLLLENLRNLEFLD FPVSKLELISQLKWVLPWLLLEEVSRSRLDIHGAVKKMQMFSEAE TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001007170
RefSeq Size:	2434
RefSeq ORF:	768
Synonyms:	ZKSCAN16; ZSCAN48
Locus ID:	158399
UniProt ID:	Q8TF39




[View online »](#)

Cytogenetics: 9q31.3

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified ZNF483 protein (Cat# [TP320978]). The protein was produced from HEK293T cells transfected with ZNF483 cDNA clone (Cat# [RC220978]) using MegaTran 2.0 (Cat# [TT210002]).