

Product datasheet for PH323130

CCDC9B (NM_207380) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	C15orf52 MS Standard C13 and N15-labeled recombinant protein (NP_997263)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC223130
Predicted MW:	57.1 kDa
Protein Sequence:	>RC223130 representing NM_207380 Red =Cloning site Green =Tags(s)


MISC AEQSRQGEAGRGAPVAPAFLLPLWLP R GCSGILSVPAVAMHSAGTPRAESPM SRQEKDAELDRRI
 VALRKKNQALLRRYQEIQEDRRQAEQGGMAVTT PALLQPDGLTVTISQVPGEKRVVSRNWARGTCGPRVT
 NEMLEDEDAEDHGGTFCLGELVELAVTMENKAEGKRIVSEKPTRARNQGIEGSPGGRVTRSPPTQVAISS
 DSARKGSWEPSRPVGPPEAGWDYAQWKQEREQIDLARLARHRDAQGDWRPWL D KAKSTLQDCSQR L
 GEGPARAGSRRGPRSHQKLQPPPLL PDGKGRGGQASRPSVAPATGSKARGKERL TGRARRWDMKEDKEEL
 EGQEGSQSTRETPSEEEQAQKQSGMEQGRLGSAPAASPALASPEGPKGESVASTASSVPCSPQEPDLAPL
 DLSLGGAGIPGPRESGCVLGLRPGAQESPVSWPEGSKQQPLGWSNHQAELEVQTCPEPQRGAGLPEPGED
 RSGKSGAQQLAPRSRPTRGGSQSRGTAGVRRRTGRPGPAGRC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_997263
RefSeq Size:	5344
RefSeq ORF:	1602


[View online »](#)

Synonyms: C15orf52
Locus ID: 388115
UniProt ID: Q6ZUT6
Cytogenetics: 15q15.1

Product images:


Coomassie blue staining of purified CCDC9B protein (Cat# [TP323130]). The protein was produced from HEK293T cells transfected with CCDC9B cDNA clone (Cat# [RC223130]) using MegaTran 2.0 (Cat# [TT210002]).