

Product datasheet for PH325433

RRAD (NM_001128850) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	RRAD MS Standard C13 and N15-labeled recombinant protein (NP_001122322)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC225433
Predicted MW:	33.1 kDa
Protein Sequence:	>RC225433 representing NM_001128850 Red=Cloning site Green=Tags(s) MTLN GGSGAGGSRGGGQERERRRGSTPWGPAPPLHRRSMPVDERDLQAALTPGAL TAAAAGTGTQGPRL DWPEDSEDSLSSGGSDSDSVYKVL LLGAPGVGKSALARIFGGVEDGPEAEAAGHTYDRSIVVDGEEASL MYYDIWEQDGGRWLPGHCMAMGDAYVIVYSVTDKGSFEKASELRVQLRRARQTDDVPIILVGNKSDLVRS REVSVDEGRACAVVFDCKFIETSAALHHNVQALFEGVVRQIRLRRDSKEANARRQAGTRRRESLGGKAKR FLGRIVARNSRKMAFRAKSKSCHDLSVL TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001122322
RefSeq ORF:	924
Synonyms:	RAD; RAD1; REM3
Locus ID:	6236
UniProt ID:	P55042 , A0A024R6X0




[View online »](#)

Cytogenetics: 16q22.1

Summary: May play an important role in cardiac antiarrhythmia via the strong suppression of voltage-gated L-type Ca(2+) currents. Regulates voltage-dependent L-type calcium channel subunit alpha-1C trafficking to the cell membrane (By similarity). Inhibits cardiac hypertrophy through the calmodulin-dependent kinase II (CaMKII) pathway. Inhibits phosphorylation and activation of CAMK2D.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified RRAD protein (Cat# [TP325433]). The protein was produced from HEK293T cells transfected with RRAD cDNA clone (Cat# [RC225433]) using MegaTran 2.0 (Cat# [TT210002]).