

Product datasheet for PH325452

Heme oxygenase 2 (HMOX2) (NM_001127205) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	HMOX2 MS Standard C13 and N15-labeled recombinant protein (NP_001120677)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC225452
Predicted MW:	36 kDa
Protein Sequence:	>RC225452 representing NM_001127205 Red=Cloning site Green=Tags(s) MSAEVETSEGVDSEKKNNGALEKENQMRMADLSELLKEGTKEAHDRAENTQFVKDFLKGNIKKELFKLA TTALYFTYSALEEEEMERNKDHPAFAPLYFPMELHRKEALTKDMEYFFGENWEEQVQCPKAAQKYVERIHY IGQNEPELLVAHAYTRYMGDLGGQVLKVAQRALKLPSTGEGTQFYLFENVDNAQQFKQLYRARMNALD LNMKTKERIVEEANKAFEYNMQIFNELDQAGSTLARETLEDGFPVHDGKGMKPCPFYAAEQDKGALEGS SCPFRTAMAVLRKPSLQFILAAGVALAAGLLAWYYM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001120677
RefSeq Size:	1776
RefSeq ORF:	948
Synonyms:	HO-2
Locus ID:	3163



[View online »](#)

UniProt ID: [P30519](#)


Cytogenetics: 16p13.3

Summary: Heme oxygenase, an essential enzyme in heme catabolism, cleaves heme to form biliverdin, which is subsequently converted to bilirubin by biliverdin reductase, and carbon monoxide, a putative neurotransmitter. Heme oxygenase activity is induced by its substrate heme and by various nonheme substances. Heme oxygenase occurs as 2 isozymes, an inducible heme oxygenase-1 and a constitutive heme oxygenase-2. HMOX1 and HMOX2 belong to the heme oxygenase family. Several alternatively spliced transcript variants encoding three different isoforms have been found for this gene. [provided by RefSeq, Oct 2013]

Protein Families: Transmembrane

Protein Pathways: Porphyrin and chlorophyll metabolism

Product images:



Coomassie blue staining of purified HMOX2 protein (Cat# [TP325452]). The protein was produced from HEK293T cells transfected with HMOX2 cDNA clone (Cat# [RC225452]) using MegaTran 2.0 (Cat# [TT210002]).