

Product datasheet for PH325515

C17orf76 (LRRC75A) (NM_001113567) Human Mass Spec Standard

Product data:


Product Type:	Mass Spec Standards
Description:	C17orf76 MS Standard C13 and N15-labeled recombinant protein (NP_001107039)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC225515
Predicted MW:	37.6 kDa
Protein Sequence:	>RC225515 representing NM_001113567 Red=Cloning site Green=Tags(s) MGRQTKGSLAERASPGAAPGPRRERPDFWASLLLRAGDKAGRAGAGMPPYHRRVGMVQELLRMVQRGR EEAGTLLQHLRQDLGMESTLDDVLYRYASFRNLVDPITHDLIIISLARYIHCPKPEGDALGAMEKLCRQL TYHLSPHSQWRRHRGLVKRKPQACLKAVLAGSPPDNTVDLSGIPLTSRDLERVTSYLQRCGEQVDSVELG FTGLTDDMVLQLPALSTLPRLTTLALNGNRLTRAVLRDLTDILKDPKSKFPNVTWIDLGNVDFISLPQP FLLSLRKRSPKQGHLPITILELGEGPGSGEEVREGTVGQEDPGGGPVAPAEDHHEGKETVAAAQT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001107039
RefSeq ORF:	1032
Synonyms:	C17orf76; FAM211A
Locus ID:	388341
UniProt ID:	Q8NAA5 , B7ZMA3



[View online »](#)

Cytogenetics: 17p11.2

Product images:



Coomassie blue staining of purified LRRC75A protein (Cat# [TP325515]). The protein was produced from HEK293T cells transfected with LRRC75A cDNA clone (Cat# [RC225515]) using MegaTran 2.0 (Cat# [TT210002]).