

Product datasheet for PH325864

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

YAP1 (NM_001130145) Human Mass Spec Standard

Product data:

Product Type: Mass Spec Standards

Description: YAP1 MS Standard C13 and N15-labeled recombinant protein (NP_001123617)

Species:HumanExpression Host:HEK293

Expression cDNA Clone

RC225864

or AA Sequence:

Predicted MW: 54.3 kDa

Protein Sequence: >RC225864 representing NM_001130145

Red=Cloning site Green=Tags(s)

MDPGQQPPPQPAPQGQGQPPSQPPQGQGPPSGPGQPAPAATQAAPQAPPAGHQIVHVRGDSETDLEALFN AVMNPKTANVPQTVPMRLRKLPDSFFKPPEPKSHSRQASTDAGTAGALTPQHVRAHSSPASLQLGAVSPG TLTPTGVVSGPAATPTAQHLRQSSFEIPDDVPLPAGWEMAKTSSGQRYFLNHIDQTTTWQDPRKAMLSQM NVTAPTSPPVQQNMMNSASGPLPDGWEQAMTQDGEIYYINHKNKTTSWLDPRLDPRFAMNQRISQSAPVK QPPPLAPQSPQGGVMGGSNSNQQQQMRLQQLQMEKERLRLKQQELLRQAMRNINPSTANSPKCQELALRS QLPTLEQDGGTQNPVSSPGMSQELRTMTTNSSDPFLNSGTYHSRDESTDSGLSMSSYSVPRTPDDFLNSV DEMDTGDTINQSTLPSQQNRFPDYLEAIPGTNVDLGTLEGDGMNIEGEELMPSLQEALSSDILNDMESVL

AATKLDKESFLTWL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Concentration: >0.05 μg/μL as determined by microplate BCA method

Labeling Method: Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3

Storage: Store at -80°C. Avoid repeated freeze-thaw cycles.

Stability: Stable for 3 months from receipt of products under proper storage and handling conditions.

RefSeq: NP 001123617

RefSeq ORF: 1512

Synonyms: COB1; YAP; YAP2; YAP65; YKI



YAP1 (NM_001130145) Human Mass Spec Standard - PH325864

Locus ID: 10413

UniProt ID: P46937, Q86T74

Cytogenetics: 11q22.1


Summary: This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is

> involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by

RefSeq, Aug 2013]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified YAP1 protein (Cat# [TP325864]). The protein was produced from HEK293T cells transfected with YAP1 cDNA clone (Cat# [RC225864]) using MegaTran 2.0

(Cat# [TT210002]).