

Product datasheet for PH325960

SH3BP2 (NM_001122681) Human Mass Spec Standard

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Mass Spec Standards |
| Description: | SH3BP2 MS Standard C13 and N15-labeled recombinant protein (NP_001116153) |
| Species: | Human |
| Expression Host: | HEK293 |
| Expression cDNA Clone or AA Sequence: | RC225960 |
| Predicted MW: | 62.2 kDa |
| Protein Sequence: | >RC225960 protein sequence Red=Cloning site Green=Tags(s) |

MAAEEMHWPVPMKAIGAQNLLTMPGGVAKAGYLHKKGGTQLQLLKWPLRFVIIHKRCVYYFKSSTSASPQ
GAFSLSGYNRMRAAEETTNNVFPFKIIHISKKHRTWFFSASSEERKSWMALLRREIGHFHEKKDLPL
DTSDDSDTDSFYGAVERPVDISLSPYPTDNEDYEHDDEDDSYLEPDSPEPGRLEDALMHPAYPPPPVP
TPRKPAFSDMPRAHSFTSKGPGPLPPPPPKHGLPDVGLAAEDSKRDPLCPRAEPCPRVPATPRRMSDP
PLSTMPTAPGLRPPCFRESASPSPEPWTGPHGACSTSSAAIMATATSRNCDKLSFHLSPRGPPTSEPP
PVPANKPKFLKIAEEDPPREAAMPGLFVPPVAPRPPALKLPVPEAMARPAVLPRPEKQPLPHLQRSPPDG
QSFRSFSFEKPRQPSQADTGGDDSDYEVKPLPNSVFNVTTESCEVERLFKATSPRGEPQDGLYICIRNS
STKSGKVLVWDETSNKVRNYRIFEKDSKFYLEGEVLFVSVGSMVEHYHTHVLPESHQSLLLRHPYGYTGP
R

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|------------------|--|
| Tag: | C-Myc/DDK |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Labeling Method: | Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3 |
| Storage: | Store at -80°C. Avoid repeated freeze-thaw cycles. |
| Stability: | Stable for 3 months from receipt of products under proper storage and handling conditions. |
| RefSeq: | NP_001116153 |
| RefSeq Size: | 9068 |
| RefSeq ORF: | 1683 |



[View online »](#)

Synonyms: 3BP-2; 3BP2; CRBM; CRPM; RES4-23

Locus ID: 6452

UniProt ID: [P78314](#), [A0A384N6E5](#)


Cytogenetics: 4p16.3

Summary: The protein encoded by this gene has an N-terminal pleckstrin homology (PH) domain, an SH3-binding proline-rich region, and a C-terminal SH2 domain. The protein binds to the SH3 domains of several proteins including the ABL1 and SYK protein tyrosine kinases, and functions as a cytoplasmic adaptor protein to positively regulate transcriptional activity in T, natural killer (NK), and basophilic cells. Mutations in this gene result in cherubism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]

Protein Families: Druggable Genome

Protein Pathways: Natural killer cell mediated cytotoxicity

Product images:



Coomassie blue staining of purified SH3BP2 protein (Cat# [TP325960]). The protein was produced from HEK293T cells transfected with SH3BP2 cDNA clone (Cat# [RC225960]) using MegaTran 2.0 (Cat# [TT210002]).