

## Product datasheet for PH326078

### SPARCL1 (NM\_001128310) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	SPARCL1 MS Standard C13 and N15-labeled recombinant protein (NP_001121782)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC226078
Predicted MW:	75.2 kDa
Protein Sequence:	>RC226078 protein sequence Red=Cloning site Green=Tags(s)

MKTGLFFLCLLGTA AAIPTNARLLSDHSPKPTAETVAPDNTAIPSLRAEDEENEKETAVSTEDDSSHKAEEK  
SSVLKSKEESHEQSAEQGKSSSQELGLKDQEDSDGDL SVNLEYAPTEGLTDIKEDMSEPQEKKLSSENTDF  
LAPGVSSFTDSNQESITKREENQEQRNYSHHQLNRSSKHSQGLRDQGNQE QDPNISNGEEEEKEPEGE  
VGT HNDNQERKTELPREHANSKQEEDNTQSDDI LEE SDQPTQVSKMQEDEFDQGNQE QEDNSNAEMEEEN  
ASNVNKHIQETEQSQEGKTGLEAISNHKETEETVSEALLMEPTDDGNTTPRNHGVDDDDGDDGGT  
DGPRHSASDDYFIPSQAFLEAERAQSIAYHLKIEEQREKVHENENIGTTEPGEHQEAKKAENSSNEEETS  
SEGNMRVHAVDSCMSFQCKRGHICKADQQGKPHVCQDPVTCPPTKPLDQVCGTDNQT YASSCHLFATKC  
RLEGTKKGHLQLDYFGACKSIPTCTDFEVIQFPLRMRDWLKNILMQLYEANSEHAGYLNEKQRNKVKKI  
YLDEKRLLAGDHPIDL LLDL RDFKKNYHMYVYPVHWQF SELDQHPMDRVLTHSELAPLRASLVPMEHCI TRF  
FEEDPNKDKHITLKEWGHCFGIKEEDIDENLLF


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001121782</a>
RefSeq Size:	3012



[View online »](#)

RefSeq ORF:	1992
Synonyms:	MAST 9; MAST9; PIG33; SC1
Locus ID:	8404
UniProt ID:	<a href="#">Q14515</a>
Cytogenetics:	4q22.1
Protein Families:	Druggable Genome, Secreted Protein

**Product images:**

Coomassie blue staining of purified SPARCL1 protein (Cat# [TP326078]). The protein was produced from HEK293T cells transfected with SPARCL1 cDNA clone (Cat# [RC226078]) using MegaTran 2.0 (Cat# [TT210002]).