

## Product datasheet for PH326168

### Breast cancer suppressor candidate 1 (VWA5A) (NM\_001130142) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	VWA5A MS Standard C13 and N15-labeled recombinant protein (NP_001123614)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC226168
Predicted MW:	86.3 kDa
Protein Sequence:	>RC226168 representing NM_001130142 Red=Cloning site Green=Tags(s)

MVHFCGLLTLHREPVPLKISISVSVNIYEFVAGVSATLNYEENEKVPLEAFFVFPMDSDAVYSFEALVDG  
KKIVAELQDKMKARTNYEKAISQGHQAFLLLEGDSSSRDVFSCNVGNLQPGSKAAVTLKYVQELPLEADGA  
LRFVLPVAVLNPRYQFSGSSKDSCLNVKTPIVPVEDLPYTL SMVATIDSQHGIEKVQSNCPSPTEYLGED  
K TSAQVSLAAGHKFDRDVELLIYNEVHTPSVLEMGMPNMPGHLMGDPSAMVSFYPNIPEDQPSNTCG  
EFIFLMDRSGSMQSPMSSQDTSQLRIQAAKETLILLKSLPIGICYFNIYGFSSYEACFPESVKYTQOTM  
EEALGRVKLMQADLGGTEILAPLQNIYRGPSIPGHPLQLFVFTDGEVTDTF SVIKEVRINRQKHRCFSFG  
IGEGTSTSLIKGIARASGGTSEFITGKDRMQSKALRTLKRSLQPVEDVSLSWHLPPGLSAKMLSPEQTV  
IFRQRLISYAQLTGRMPAAETTGEVCLKYTLQGKTFEDKVTFPLQPKPDVNLTIHRLAAKSLLTQTKDMG  
LRETPASDKKDALNLSLESGVISSFTAFIAINKELNKPVQGPLAHRDVRPILLGASAPLKIKCQSGFRK  
ALHSDRPPSASQPRGELMCYKAKTFQMDDYSLCGLISHKDQHSFGFGENHLVQLIYHQANGSWDLNEDL  
AKILGMSLEEIMAAQPAELVDSSGWATILAVIWLHNSGKDLKCEWELLERKAVAWMRAHAGSTMPSSVKA  
AITFLKSSVDP AIFAF

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>4</sub> ]-L-Arginine and [U- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> ]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001123614</a>



[View online »](#)

RefSeq ORF:	2358
Synonyms:	BCSC-1; BCSC1; LOH11CR2A
Locus ID:	4013
UniProt ID:	<a href="#">O00534</a> , <a href="#">A0A024R3H3</a>
Cytogenetics:	11q24.2
Summary:	May play a role in tumorigenesis as a tumor suppressor. Altered expression of this protein and disruption of the molecular pathway it is involved in, may contribute directly to or modify tumorigenesis.[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified VWA5A protein (Cat# [TP326168]). The protein was produced from HEK293T cells transfected with VWA5A cDNA clone (Cat# [RC226168]) using MegaTran 2.0 (Cat# [TT210002]).