

Product datasheet for PH326773


FAM76A (NM_001143914) Human Mass Spec Standard

Product data:

Product Type:	Mass Spec Standards
Description:	FAM76A MS Standard C13 and N15-labeled recombinant protein (NP_001137386)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC226773
Predicted MW:	31.7 kDa
Protein Sequence:	>RC226773 representing NM_001143914 Red=Cloning site Green=Tags(s) MAALYACTKCHQRFPEALSQGGQLCKECRIAHVPVKCTYCRTEYQQESKTNTICKKCAQNVQLYGTPKP CQYCNIIAAFIGNKCQRCTNSEKKYPPYSCEQCKQQCAFDKDDRKVDGKLLCWLCTLSYKRVLQKTK EQRKHLSSSRAGHQEKEQYSRLSGGGHYNSFSPDLALDSPGTDHFVIIAQLKEEVATLKKMLHQKDQMI LEKEKKITELKADFQYQESQMRAMNQMEKTHKEVTEQLQAKNRELLKQAAALSKSKKSEKSGAITSP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- ¹³ C ₆ , ¹⁵ N ₄]-L-Arginine and [U- ¹³ C ₆ , ¹⁵ N ₂]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	NP_001137386
RefSeq ORF:	834
Locus ID:	199870
UniProt ID:	Q8TAV0
Cytogenetics:	1p35.3


[View online »](#)

Product images:



Coomassie blue staining of purified FAM76A protein (Cat# [TP326773]). The protein was produced from HEK293T cells transfected with FAM76A cDNA clone (Cat# [RC226773]) using MegaTran 2.0 (Cat# [TT210002]).