

## Product datasheet for PH326856

### TMEM51 (NM\_001136218) Human Mass Spec Standard

#### Product data:

Product Type:	Mass Spec Standards
Description:	TMEM51 MS Standard C13 and N15-labeled recombinant protein (NP_001129690)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	RC226856
Predicted MW:	27.8 kDa
Protein Sequence:	>RC226856 protein sequence Red=Cloning site Green=Tags(s)  MMAQSKANGSHYALTAIGLGMLVLGVIMAMWNLVPGFSAAEKPTAQGSNKTEVGGGILKSKTFSVAYVLV GAGVMLLLLSICLSIRDKRKQRQGEDLAHVQHPTGAGPHAQEEDSQEEEEDEEAASRYVPSYEEVMNT NYSEARGEEQNPRLSISLPSYESLTGLDETTPTSTRADVEASPGNPPDRQNSKLAKRLKPLKVRRIKSEK LHLKDFRINLPDKNVPPPSIEPLTPPPQYDEVQEKAPDTRPPD  TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Labeling Method:	Labeled with [U- 13C6, 15N4]-L-Arginine and [U- 13C6, 15N2]-L-Lysine
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3
Storage:	Store at -80°C. Avoid repeated freeze-thaw cycles.
Stability:	Stable for 3 months from receipt of products under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001129690</a>
RefSeq Size:	2000
RefSeq ORF:	759
Synonyms:	C1orf72
Locus ID:	55092
UniProt ID:	<a href="#">Q9NW97</a> , <a href="#">A0A024QZ97</a>




[View online »](#)

Cytogenetics: 1p36.21

Protein Families: Transmembrane

### Product images:



Coomassie blue staining of purified TMEM51 protein (Cat# [TP326856]). The protein was produced from HEK293T cells transfected with TMEM51 cDNA clone (Cat# [RC226856]) using MegaTran 2.0 (Cat# [TT210002]).