

Product datasheet for **PP021B2**

6Ckine / CCL21a Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA, WB

Recommended Dilution: **Elisa**

Direct: To detect Mouse Exodus-2 by direct ELISA (using 100 µl/well antibody solution) a concentration of 0.25–1.0 µg/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with compatible secondary reagents, allows the detection of at least 0.2–0.4 ng/well of recombinant Mouse Exodus-2.

Sandwich: To detect Mouse Exodus-2 by sandwich ELISA (using 100 µl/well antibody solution) a concentration of 0.25–1.0 µg/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with Anti-Murine Exodus-2 (PP021P1 or PP021P2) as a capture antibody, allows the detection of at least 0.2–0.4 ng/well of recombinant Mouse Exodus-2.

Western Blot: To detect mExodus-2 by Western Blot analysis this antibody can be used at a concentration of 0.1–0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant Mouse Exodus-2 is 1.5–3.0 ng/lane, under either reducing or non-reducing conditions.

Reactivity: Mouse

Host: Rabbit

Clonality: Polyclonal

Immunogen: Highly pure (>98%) recombinant Mouse Exodus-2.

Specificity: This antibody recognises Murine Exodus-2.

Formulation: PBS, pH 7.2 without preservatives.

Label: Biotin

State: Lyophilized (sterile filtered) purified Ig fraction.

Reconstitution Method: Restore in sterile PBS containing 0.1% BSA to a concentration of 0.1–1.0 mg/ml.

Purification: Affinity Chromatography.

Conjugation: Biotin

Storage: Store the antibody prior to reconstitution at -20°C. Following reconstitution the antibody can be stored at 2–8°C for one month or at -20°C for longer.
Avoid repeated freezing and thawing.



[View online »](#)

Stability:	Shelf life: One year from despatch.
Gene Name:	C-C motif chemokine 21c
Database Link:	Entrez Gene 100041504 Mouse P86792
Background:	<p>Exodus 2 is one of several CC cytokine genes clustered on the p arm of chromosome 9. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterised by two adjacent cysteines. Similar to other chemokines the protein encoded by this gene inhibits hemopoiesis and stimulates chemotaxis. This protein is chemotactic in vitro for thymocytes and activated T cells, but not for B cells, macrophages, or neutrophils. The cytokine encoded by this gene may also play a role in mediating homing of lymphocytes to secondary lymphoid organs. It is a high affinity functional ligand for chemokine receptor 7 (CCR7) that is expressed on T and B lymphocytes and a known receptor for another member of the cytokine family (small inducible cytokine A19).</p>
Synonyms:	C-C motif chemokine 21a, Beta chemokine exodus-2, Small-inducible cytokine A21a, Scya21a, TCA4, TCA-4
Note:	Centrifuge vial prior to opening!