

# Product datasheet for PP1041B2

## **Leptin (LEP) Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type: Primary Antibodies** 

**Applications:** ELISA, WB

Recommended Dilution: Direct ELISA: To detect Human Leptin by Direct ELISA (using 100 µl/well antibody solution) a

concentration of 0.25-1.0 µg/ml of this antibody is required. This antibody, in conjunction with

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

compatible secondary reagents, allows the detection of at least 0.2-0.4 ng/well of

recombinant Human Leptin.

Sandwich ELISA: To detect Human Leptin by Sandwich ELISA (using 100 µl/well antibody solution) a concentration of 0.25-1.0 µg/ml of this antibody is required. This biotinylated antibody, in conjunction with Anti-Human Leptin (Cat.-No PP1041P) as a capture antibody,

allows the detection of at least 0.2-0.4 ng/well of recombinant hLeptin.

Western Blot: To detect Human Leptin by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant Human Leptin is 1.5-3.0 ng/lane, under either reducing or

non-reducing conditions.

Reactivity: Human Host: Rabbit

Polyclonal Clonality:

Highly pure (>98%) E.coli derived recombinant Human Leptin (Cat.-No PA094) Immunogen:

Specificity: This antibody reacts with Human Leptin.

Other species not tested.

Formulation: PBS, pH 7.2 without preservatives

Label: Biotin

State: Lyophilized (sterile filtered) purified Ig fraction.

Label: conjugated

**Reconstitution Method:** Restore in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.

Centrifuge vial prior to opening!

**Purification:** Affinity Chromatography

Conjugation: Biotin





#### Leptin (LEP) Rabbit Polyclonal Antibody - PP1041B2

**Storage:** Prior to reconstitution store at 2-8°C.

Following reconstitution store undiluted at 2-8°C for one month

or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

Gene Name: leptin

**Database Link:** Entrez Gene 3952 Human

P41159

**Background:** Leptin plays a critical role in the regulation of body weight by inhibiting food intake and


stimulating energy expenditure. Defects in Leptin production cause severe hereditary obesity in rodents and humans. In addition to its effects on body weight, leptin has a variety of other functions, including the regulation of hematopoiesis, angiogenesis, wound healing, and the immune and inflammatory response. The Leptin gene is the human homolog of the gene (ob)

mutant in the mouse 'obese' phenotype.


Defects in the Leptin gene are the cause of profound obesity and type II diabetes.

Synonyms: LEP, OB, OBS, Obesity factor, Obese protein

### **Product images:**



Direct ELISA using Biotin conjuagted Leptin Antibody Cat.-No [PP1041B]



Sandwich ELISA using Biotin conjuagted Leptin Antibody Cat.-No [PP1041B]