

Product datasheet for **PP1203P2**

BAFF (TNFSF13B) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Neutralization. ELISA: To detect hBAFF by direct ELISA (using 100 µl/well antibody solution) a concentration of at least 0.5 µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2-0.4 ng/well of recombinant hBAFF. Western Blot: To detect hBAFF by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hBAFF is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Highly pure (>98%) recombinant hBAFF.
Formulation:	0.5 X PBS, pH 7.4 without preservatives. State: Aff - Purified State: Lyophilized purified Ig fraction.
Reconstitution Method:	Restore in sterile water to a concentration of > 0.2 mg/ml.
Purification:	Affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody prior to reconstitution at -20°C. Following reconstitution the antibody can be stored at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	tumor necrosis factor superfamily member 13b
Database Link:	Entrez Gene 10673 Human Q9Y275



[View online »](#)

Background:

Members of the TNF superfamily regulate immune responses and induce apoptosis. A novel member in the TNF family was recently identified by several groups and designated BAFF (for B cell Activating Factor belonging to the TNF Family), BLYS (for B Lymphocyte Stimulator), TALL1 (for TNF- and ApoL- related Leukocyte-expressed Ligand), and THANK (for TNF Homologue that Activate Apoptosis, NFkB and c-jun N-terminal Kinase). BAFF/BLYS was characterized as a B cell stimulator since it induced B cell proliferation and immunoglobulin secretion. Two receptors for BAFF were recently identified and designated TACI and BCMA. BAFF also signals through a third TNF receptor BAFFR/BR3. BAFF and its receptors are involved in the development of systemic lupus erythaematosus and other B cell associated autoimmune diseases. Like TNFa and TRAIL, THANK was shown to activate NF-kB and c-jun N terminal kinase (JNK) and to induce apoptosis.

The human BAFF gene codes for a 285 amino acid type II transmembrane protein containing a 46 amino acid cytoplasmic domain, a 21 amino acid transmembrane domain, and a 218 amino acid extracellular domain.

Synonyms:

TNFSF13B, BLYS, TALL1, TNFSF20, ZTNF4