

Product datasheet for **PP1212P1**

LIX (CXCL6) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Sandwich ELISA: To detect mLIX by sandwich ELISA (using 100 µl/well antibody solution) a concentration of 0.5-2.0 µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Murine LIX (PP1212B) as a detection antibody, allows the detection of at least 0.2-0.4 ng/well of recombinant mLIX. Western Blot: To detect mLIX, this antibody can be used at a concentration of 0.1-0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant mLIX is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Highly pure (>98%) <i>E.coli</i> derived Recombinant Murine LIX (CXCL6) (92 a.a.).
Formulation:	0.5 X PBS pH 7.4 without preservatives. State: Aff - Purified State: Lyophilized purified Ig fraction.
Reconstitution Method:	Restore in sterile water to a concentration of > 0.2 mg/ml.
Purification:	Affinity chromatography employing immobilized mLIX matrix.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	C-X-C motif chemokine 6, Small-inducible cytokine B6, Granulocyte chemotactic protein 2, Chemokine alpha 3, SCYB6, CKA-3, CKA-3



[View online »](#)

Product images:

