

## Product datasheet for **PP1221P2**

### Adiponectin (ADIPOQ) (N-term Variant) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: To detect hAdipolean Variant by direct ELISA (using 100 µl/well antibody solution) a concentration of at least 0.5µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2 - 0.4 ng/well of recombinant hgAcrp30/Adipolean. Western Blot: To detect hAdipolean Variant by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hAdipolean is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Highly purified (>98%) recombinant hAdipolean Variant.
Formulation:	0.5 x PBS, pH 7.4 without preservatives. State: Aff - Purified State: Lyophilized purified Ig fraction.
Reconstitution Method:	Restore in sterile water to a concentration of > 0.2 mg/ml.
Purification:	Affinity chromatography employing immobilized hAdipolean Variant matrix.
Conjugation:	Unconjugated
Storage:	Store the antibody prior to reconstitution at -20°C. Following reconstitution the antibody can be stored at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	adiponectin, C1Q and collagen domain containing
Database Link:	<a href="#">Entrez Gene 9370 Human Q15848</a>
Synonyms:	ADIPOQ, ACDC, ACRP30, APM1, GBP28



[View online »](#)