

Product datasheet for **PP1225B2**

RELM beta (RETNLB) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: To detect hRELM β by direct ELISA (using 100 μ l/well antibody solution) this antibody can be used at a concentration of 0.15 - 0.30 μ g/ml. Used in conjunction with compatible secondary reagents, allows the detection of at least 0.2 ng/well of recombinant hRELM β . Western blot: To detect hRELM β by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 μ g/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hRELM β is 1.5 -3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Highly pure (>98%) recombinant hRELM β (human RELM-beta).
Formulation:	BPS, pH 7.2 without preservatives. Label: Biotin State: Lyophilized purified fraction. Label: Conjugated to
Purification:	Affinity chromatography.
Conjugation:	Biotin
Storage:	Store the antibody prior to reconstitution at -20°C. Following reconstitution the antibody can be stored at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	resistin like beta
Database Link:	Entrez Gene 84666 Human Q9BQ08
Synonyms:	Resistin-like beta, RELMbeta, RELM-beta RETNLB, CCRG, FIZZ2, HXCP2, RETNL2



[View online »](#)