

Product datasheet for **PP1225P1**

RELM beta (RETNLB) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Immunostaining. ELISA: To detect hRELM β by direct ELISA (using 100 μ l/well antibody solution) a concentration of at least 0.5 μ g/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2-0.4 ng/well of recombinant hRELM β . Western blot: To detect hRELM β by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 μ g/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hRELM β is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Highly pure (>98%), E.coli derived, 19.0 kDa recombinant Human RELM β (human RELM-beta).
Formulation:	BPS, pH 7.2 without preservatives. State: Aff - Purified State: Lyophilized (sterile filtered) purified Ig fraction
Reconstitution Method:	Restore in sterile water to a concentration of 0.1-1.0 mg/ml.
Purification:	Affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	resistin like beta
Database Link:	Entrez Gene 84666 Human Q9BQ08



[View online »](#)

Synonyms:

Resistin-like beta, RELMbeta, RELM-beta RETNLB, CCRG, FIZZ2, HXCP2, RETNL2