

Product datasheet for **PP1228P1**

TL1A (TNFSF15) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Neutralisation: Data not available. ELISA: To detect hTL-1A by direct ELISA (using 100 µl/well antibody solution) a concentration of at least 0.5 -2.0 µg/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2 - 0.4 ng/well of recombinant hTL-1A. Western blot: To detect hTL-1A by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hTL-1A is 1.5 -3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Highly pure (>98%) E.coli derived 22.0 kDa recombinant hTL-1A (human TL-1A).
Formulation:	PBS, 7.4 without preservatives. State: Aff - Purified State: Lyophilized purified fraction.
Reconstitution Method:	Restore in sterile water to a concentration of 0.1-1.0 mg/ml.
Purification:	Affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	tumor necrosis factor superfamily member 15
Database Link:	Entrez Gene 9966 Human O95150



[View online »](#)

Synonyms: TL1, TL1A, VEGI