



Product datasheet for PR100001

pRMT-Luc, Promoter-less Firefly Luciferase Reporter Vector

Product data:

- Product Type:** Reporter Vectors
- E. coli Selection:** Ampicillin
- Features:**
- Versatile MCS for promoter sequence cloning
 - Low background
 - Neo selection marker
- Reporter Gene:** Firefly Luciferase

Schematic of the multiple cloning sites:




Storage: Store the vector DNA at -20°C.

Stability: Vector DNA is stable for at least one year from date of receipt when stored at -20°C.



[View online »](#)

Product images:



pRMT-Luc

```

Nhe I          EcoRI          BamHI  Kpn I
GCTAGCTCTAGAAGAAATTTTGTAAACGACTCACTATAGGGCGCGGGGAATTCGCTGACTGGATCCGGTACC

Sgf I  Asc I  Hind III  Spe I  Mlu I
GAGGAGATCTGCCCGCGATCGCCGCGCGCCAGATCTCAAGCTTAACTAGTTAGCGGACCGACGGCTAC

Not I  Xho I
GCGGCCGCTCGAG
    
```

