


Product datasheet for **PR100002**

pRMT-tGFP, Promoter-less GFP Reporter Vector

Product data:

- Product Type:** Reporter Vectors
- E. coli Selection:** Ampicillin (100 ug/mL)
- Features:**
- Versatile MCS for promoter sequence cloning
 - Low background
 - Neo selection marker
- Reporter Gene:** tGFP

Schematic of the multiple cloning sites:




- Storage:** Store the vector DNA at -20°C.
- Stability:** Vector DNA is stable for at least one year from date of receipt when stored at -20°C.



[View online »](#)

Product images:



pRMT-tGFP

Nhe I EcoRI BamH I Kpn I
GCTAGCTCTAGAAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACC

Sgf I Asc I Hind III Spe I Mlu I
GAGGAGATCTGCCCGCGATCGCCGGCGGCCAGATCTCAAGCTTAAGTTAGTCGGACCGACGCGTAC

Not I Xho I
GCGGCCGCTCGAG