

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for PR300004

5x Passive lysis buffer

Product data:

Product Type:	Luciferase Assay Kits
Description:	5x Passive lysis buffer
Components:	1 bottle of 30 ml Passive lysis buffer for use with Renilla Luciferase Assay Kits
Features:	Concentrated lysis buffer designed for use with our Renilla Luciferase Assays
Assay Type:	Accessory
Summary:	5x Passive Lysis Buffer is concentrated lysis buffer designed for use with Renilla luciferase assays (Cat# PR300002 & PR300007).
Storage:	Store at -20°C or below.
Stability:	The buffer is stable for at least 12 months from date of receipt when stored at -20°C or below.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2021 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US