

Product datasheet for **PR300004**

5x Passive lysis buffer

Product data:

Product Type:	Luciferase Assay Kits
Description:	5x Passive lysis buffer
Components:	1 bottle of 30 ml Passive lysis buffer for use with Renilla Luciferase Assay Kits
Features:	Concentrated lysis buffer designed for use with our Renilla Luciferase Assays
Assay Type:	Accessory
Summary:	5x Passive Lysis Buffer is concentrated lysis buffer designed for use with Renilla luciferase assays (Cat# PR300002 & PR300007).
Storage:	Store at -20°C or below.
Stability:	The buffer is stable for at least 12 months from date of receipt when stored at -20°C or below.



[View online »](#)