

Product datasheet for **PR300005**

5x Firefly luciferase lysis buffer

Product data:

Product Type:	Luciferase Assay Kits
Description:	5x Firefly Luciferase lysis buffer
Components:	1 bottle of 15 ml Firefly luciferase lysis buffer for use with Firefly Luciferase Assay Kits
Features:	Concentrated lysis buffer designed for use with our Firefly Luciferase Assays
Assay Type:	Accessory
Summary:	5x Firefly Luciferase lysis buffer is concentrated lysis buffer designed for use with Firefly luciferase assay (Cat# PR300001 and PR300006).
Note:	This Buffer is not compatible with Renilla Luciferase. If you require additional lysis buffer for use with Renilla Luciferase Assay Kit (Cat# PR300002 and PR300007), order Passive Lysis Buffer (Cat# PR300004).
Storage:	Store at -20°C or below.
Stability:	The buffer is stable for at least 12 months from date of receipt when stored at -20°C or below.



[View online »](#)