

## Product datasheet for **PR300010**

### 1x Passive lysis buffer 2.0

#### Product data:

<b>Product Type:</b>	Luciferase Assay Kits
<b>Description:</b>	1x Passive lysis buffer 2.0
<b>Components:</b>	1 bottle of 100 ml 1x Passive lysis buffer 2.0
<b>Features:</b>	For use with Firefly & Renilla Luciferase Single Tube Assay Kit (Cat # PR300008 and PR3000009)
<b>Assay Type:</b>	Accessory
<b>Summary:</b>	1X Passive Lysis Buffer 2.0 is a cell lysis buffer that is designed for use with Firefly & Renilla Luciferase Single Tube Assay Kit (PR300008 and PR300009). Those kits are supplied with sufficient lysis buffer to perform the stated number of assays with cells grown in 96 – 24 well plates. 1X Passive Lysis Buffer 2.0 is sold separately for use with the kits when more lysis buffer (e.g. >100 uL/well) is required. The buffer also can be used to dilute recombinant firefly or Renilla luciferase enzymes without the need to add BSA to stabilize the enzymes.
<b>Storage:</b>	Store at -20°C or below.
<b>Stability:</b>	The buffer is stable for at least 12 months from date of receipt when stored at -20°C or below.



[View online »](#)