

## Product datasheet for **PR300010**

### 1x Passive lysis buffer 2.0

#### Product data:

Product Type:	Luciferase Assay Kits
Description:	1x Passive lysis buffer 2.0
Components:	1 bottle of 100 ml 1x Passive lysis buffer 2.0
Features:	For use with Firefly & Renilla Luciferase Single Tube Assay Kit (Cat # PR300008 and PR3000009)
Assay Type:	Accessory
Summary:	1X Passive Lysis Buffer 2.0 is a cell lysis buffer that is designed for use with Firefly & Renilla Luciferase Single Tube Assay Kit (PR300008 and PR300009). Those kits are supplied with sufficient lysis buffer to perform the stated number of assays with cells grown in 96 – 24 well plates. 1X Passive Lysis Buffer 2.0 is sold separately for use with the kits when more lysis buffer (e.g. >100 uL/well) is required. The buffer also can be used to dilute recombinant firefly or Renilla luciferase enzymes without the need to add BSA to stabilize the enzymes.
Storage:	Store at -20°C or below.
Stability:	The buffer is stable for at least 12 months from date of receipt when stored at -20°C or below.



[View online »](#)