

Product datasheet for PS100133

pCMV6-AC-tGFP-rtTA-Hyg Inducible Expression Vector

Product data:

Product Type: Inducible Expression Vectors

C-Tag: tGFP


E. coli Selection: Ampicillin (100 ug/mL)

Mammalian Cell Selection: Hygromycin

Features: This vector is a tetracycline-inducible mammalian expression vector. ORFs cloned in this vector will be expressed under tetracycline-inducible promoter as a fusion protein with the C-terminal tGFP tags.

- TRE_F Primer (GE100141) can be used to sequence the target sequence after cloning.
- tGFP_R Primer (GE100142) can be used to sequence the target sequence after cloning.

Product images:



Plasmid Map of pCMV6-AC-tGFP-rtTA-Hyg Inducible Expression Vector.



[View online »](#)