

## Product datasheet for **QP100003**

### **SYBR Master Mix Hi-ROX (2x)**

#### **Product data:**

<b>Product Type:</b>	qPCR Reagents
<b>Components:</b>	Sixteen vials of QP100001S each containing 1.25 mL SensiMix SYBR Master Mix - total sufficient for 800 rxns
<b>Quality Control:</b>	The SensiMix SYBR Hi-ROX Kit has been developed for highly accurate real-time PCR and has been validated on real-time PCR platforms that require high concentrations of the passive reference dye ROX for normalization of well-to-well differences.
<b>Storage:</b>	Upon receipt, please store the Master Mix at -20°C.
<b>Stability:</b>	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)