

Product datasheet for **QP100004**

SYBR Master Mix Low-ROX (2x)

Product data:

Product Type:	qPCR Reagents
Components:	One vial of QP100004S containing 1.25 mL SensiMix SYBR Master Mix – Low Rox Kit (sufficient for 100 rxns) One Datasheet
Quality Control:	The SensiMix SYBR Low-ROX Kit has been developed for highly accurate real-time PCR and has been validated on real-time PCR platforms that require low concentrations of the passive reference dye ROX for normalization of well-to-well differences.
Storage:	Upon receipt, please store the Master Mix at -20°C.
Stability:	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)