

## Product datasheet for **QP100006**

### SYBR Master Mix Low-ROX

#### Product data:

|                         |  |
|-------------------------|--|
| <b>Product Type:</b>    | qPCR Reagents  |
| <b>Quality Control:</b> | The SensiMix SYBR Low-ROX Kit has been developed for highly accurate real-time PCR and has been validated on real-time PCR platforms that require low concentrations of the passive reference dye ROX for normalization of well-to-well differences. |
| <b>Storage:</b>         | Upon receipt, please store the Master Mix at -20°C.  |
| <b>Stability:</b>       | The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.  |



[View online »](#)