

Product datasheet for **QP100009**

SYBR Master Mix Fluorescein (2x)

Product data:

Product Type:	qPCR Reagents
Components:	Sixteen vials of QP100007S each containing 1.25 mL SensiMix SYBR Master Mix – Fluorescein Kit (total sufficient for 1600 rxns) One Datasheet
Quality Control:	The SensiMix SYBR & Fluorescein Kit has been developed for highly accurate real-time PCR and has been validated on real-time PCR platforms that require the passive reference dye fluorescein for collection of dynamic well factors.
Storage:	Upon receipt, please store the Master Mix at -20°C.
Stability:	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)