

## Product datasheet for **QP100010**

### SYBR Master Mix No-Rox (2x)

#### Product data:

<b>Product Type:</b>	qPCR Reagents
<b>Components:</b>	One vial of QP100010S containing 1.25 mL SensiMix SYBR Master Mix - No-ROX Kit (sufficient for 100 rxns) One Datasheet
<b>Quality Control:</b>	The SensiMix SYBR No-ROX Kit has been developed for highly accurate real-time PCR and has been validated on real-time PCR platforms that do not require a passive reference dye for normalization of well-to-well differences.
<b>Storage:</b>	Upon receipt, please store the Master Mix at -20°C.
<b>Stability:</b>	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)