

Product datasheet for **QP100012**

SYBR Master Mix No-Rox (2x)

Product data:

Product Type:	qPCR Reagents
Components:	Sixteen vials of QP100010S each containing 1.25 mL SensiMix SYBR Master Mix (total sufficient for 1600 rxns) One Datasheet
Quality Control:	The SensiMix SYBR No-ROX Kit has been developed for highly accurate real-time PCR and has been validated on real-time PCR platforms that do not require a passive reference dye for normalization of well-to-well differences.
Storage:	Upon receipt, please store the Master Mix at -20°C.
Stability:	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)