

Product datasheet for **QP100015**

SensiFAST SYBR Master Mix Hi-ROX (2x)

Product data:

Product Type:	qPCR Reagents
Components:	Ten vials of QP100013S each containing 1.0 mL SensiFAST SYBR Master Mix - Hi-ROX Kit (sufficient for 1000 rxns) One Datasheet
Quality Control:	The SensiFAST SYBR Hi-ROX Kit has been developed for fast, highly accurate real-time PCR and has been validated on all commonly-used real-time instruments that require a high concentration of the passive reference dye ROX. A combination of the latest advances in buffer chemistry and PCR enhancers ensures that the SensiFAST SYBR Hi-ROX Kit produces reliable assay results under fast thermal cycling conditions. An antibody-mediated hot-start DNA polymerase system promotes highly-specific amplification, in turn improving assay sensitivity and dynamic range.
Storage:	Upon receipt, please store the Master Mix at -20°C.
Stability:	The Master Mix is stable for six months from date of shipping when stored at -20°C. Avoid repeated freeze-thaws as this will degrade the product.



[View online »](#)